



How to reshape the PC101 & PC102 Ergonomic U-Shaped Joystick Handle to better fit the user's hand for a more secure grip



Standard Ergonomic U-Shaped Joystick Handle



Modified Ergonomic U-Shaped Joystick Handle
Right and left sides bent inward

Immerse end of handle in boiling water (100° C) for 3 minutes. Do not immerse stem or center logo area. Remove, bend to shape, and hold until cool.



Repeat for second side, if needed.

This technique is effective for moderate bends, and will not damage the material, which can be heated and reheated to modify the shape. For more significant changes to the shape consult with an orthotist.