



Changing
the **Shape** of
Support

 Bodypoint[®]

Introducing the next advancement in shoulder harnesses

The **PivotFit™ Shoulder Harness** represents the latest in innovative design from Bodypoint. **Distributing pressure across a greater surface area** than our traditional H-style Harness, the PivotFit Shoulder Harness offers a body-contouring shape for **optimal support** and a **comfortable fit**.




Patented dual-zone elasticity distributes stretch and pressure evenly for increased comfort and support (PivotFit™ Dynamic only).

Unique shape contours to the shoulders and ribcage for a more supportive, comfortable fit.

New, one-of-a-kind, high-strength swivel buckle allows the pads to pivot into the best placement for the user's body shape and comfort.

Measuring for Bodypoint Harnesses



Measure shoulder width. If the shoulder width falls between two sizes, consider other factors such as growth, weight changes and clothing.

XL	48–56cm (19–22")
L	41–48cm (16–19")
M	33–41cm (13–16")
S	28–33cm (11–13")
XS	24–28cm (9-1/2–11")
PivotFit Harness only:	
XXS	21–24cm (8-1/4–9-1/2")

Available in sizes for 2-4 year olds

SKU's

Standard PivotFit™ Shoulder Harness

Front-Pull

SH250XL	Extra large
SH250L	Large
SH250M	Medium
SH250S	Small
SH250XS	Extra small
SH250XXS	Extra extra small

Rear-Pull

SH260XL	Extra large
SH260L	Large
SH260M	Medium
SH260S	Small
SH260XS	Extra small
SH260XXS	Extra extra small

PivotFit™ Dynamic Shoulder Harness


Front-Pull

SH270XL	Extra large
SH270L	Large
SH270M	Medium
SH270S	Small
SH270XS	Extra small
SH270XXS	Extra extra small

Rear-Pull


SH280XL	Extra large
SH280L	Large
SH280M	Medium
SH280S	Small
SH280XS	Extra small
SH280XXS	Extra extra small

Versatile attachment options for maximum adjustability




Front-pull attachment

Offering maximum stability, the front-pull attachment allows caregivers to adjust the harness while facing the user or from behind.



Rear-pull attachment

The rear-pull attachment enables caregivers to adjust the harness from behind and keeps the buckle away from the user. It's also ideal for thicker backrests.



Rear-pull option

When cam buckles are preferred for frequent length adjustment, the side release buckles can be removed from the top straps for this simple conversion.

Mounting Options

Complete your product sku by choosing -B1 or -B2 (ex. SH250XL-B1)

Flat-Mount™



Cinch-Mount™



Cam buckle also available

Sold separately.



Sold as part number FS032-2 (pair of 1" cam buckles, suitable for any model of harness)



STRENGTHEN
YOUR
POSITION

Learn more today at www.bodypoint.com/PivotFit

1.206.405.4555 | Seattle, USA | © 2013 Bodypoint. All Rights Reserved. BMM059 April 2013 Rel 1