

A New Concept in Pelvic Control


SUB-ASIS PADS AND BELTS

Molded SynerGel™-filled pad construction offers firm, yet comfortable support




SUB-ASIS PADS combined with a special hook-and-loop-compatible belt, provide two firm points of contact with the pelvis. The result is greatly enhanced stability and improved control of obliquity and rotation. The unique SynerGel™-filled pad design comfortably spaces the belt away from the abdomen, reducing pressure on the bladder.

SELECTING HIP BELTS: Measure hip width across the greater trochanters with the person seated. If the hip width falls between sizes on our chart, consider other factors such as growth, weight changes, and clothing.



XS	5"-9" (13cm-23cm)
S	7"-11" (18cm-28cm)
M	9"-15" (23cm-38cm)
L	13"-19" (33cm-48cm)



Hook and loop attachment allows precise Sub-ASIS pad placement.



Shaped pads fit the pelvic Sub-ASIS notch, providing stability while minimizing abdominal compression.

CENTER PULL | TWO-POINT

Rehab Latch™

LB225M-A2	Med, Framesaver
LB225M-B2	Med, w/o Framesaver
LB225L-A2	Lg, Framesaver
LB225L-B2	Lg, w/o Framesaver

Push Button

LB205M-A2	Med, Framesaver
LB205M-B2	Med, w/o Framesaver

REAR PULL | TWO-POINT

Rehab Latch™

LB229M-A2	Med, Framesaver
LB229M-B2	Med, w/o Framesaver

CENTER-PULL | MEDIAL FOUR-POINT

Push Button

LB405M-A2	Med, Framesaver
LB405M-B2	Med, w/o Framesaver
LB405L-A2	Lg, Framesaver
LB405L-B2	Lg, w/o Framesaver

CENTER-PULL | MEDIAL FOUR-POINT

Rehab Latch™

LB425M-A2	Med, Framesaver
LB425M-B2	Med, w/o Framesaver
LB425L-A2	Lg, Framesaver
LB425L-B2	Lg, w/o Framesaver

Extra Small Metal Push-Button

LB435XS-A2	XS, Framesaver
LB435XS-B2	XS, w/o Framesaver

DUAL-PULL | MEDIAL FOUR-POINT

Side Release

LB416XS-A2	XS, Framesaver
LB416XS-B2	XS, w/o Framesaver

REAR-PULL | MEDIAL FOUR-POINT

Rehab Latch™

LB429M-A2	Med, Framesaver
LB429M-B2	Med, w/o Framesaver

SUB-ASIS PADS

SynerGel, Velcro compatible, Pkg. of 2

PD402-2	3 3/4" (95mm)
PD401-2	4" (102mm)

