

NEW!

Joystick Handle Display and Evaluation Tools

There are now two options for showcasing Bodypoint® joystick handles. Between our updated seven-pin Joystick Handle Tabletop Display and our all-new, portable Joystick Handle Ring, your joystick handle evaluation and display needs are covered.



BMR-27



BMR-45

Joystick Handle Tabletop Display

- Durable, stationary display
- Perfect for use in clinics, hospitals, showrooms, and trade shows

NEW! Joystick Handle Ring

- Lightweight and portable
- Perfect for evaluations, trainings, and other activities in the field
- Tags include QR code linked to Bodypoint Joystick Handle webpage

Shared Features




Kit Contents

Includes one each of all Bodypoint joystick handles, plus two extra Grooved Mushroom Joystick Handles (PC110) cut down to different heights.

Removable Handles

Joystick handles easily pull off from metal pins. Convenient for demonstrations, fittings and evaluations.



Product Information

A branded tag on the Joystick Handle Ring and a label on the Joystick Tabletop Display include part numbers and product names.

Ø1/4" stem size is used on both displays.